

# CA<sup>®</sup> TIP Pliers (CA<sup>®</sup> TIP 1-6)

## Description:

### CA<sup>®</sup> TIP 1

Pliers for aligners with large spherical tip allow insertion of retention points in CA<sup>®</sup> splints, for increase of friction.

### CA<sup>®</sup> TIP 2

Pliers for aligners with small spherical tip allow forming of activation points in CA<sup>®</sup> splints.

### CA<sup>®</sup> TIP 3

Pliers for aligners with oval plateau, enable creation of activation surfaces in CA<sup>®</sup> splints.

### CA<sup>®</sup> TIP 4

Pliers for aligners destined for torque control with curved plateau, enable creation of labial / lingual activation surfaces in CA<sup>®</sup> splints.

### CA<sup>®</sup> TIP 5

Pliers for aligners with angled pliers head for correction of splint edges in case of gingival pressure marks.

### CA<sup>®</sup> TIP 6

Pliers for aligners for removal of activation points and surfaces inserted by means of CA<sup>®</sup> TIP 2, 3 and 4. Also suitable for finishing and smoothing of the splint surface.

## Cleaning/Care:

We recommend thermal disinfection (delicate program not exceeding 121°C) or spray disinfection.

## Warranty:

5 years on material, workmanship, breakage and corrosion.  
All pliers are made of stainless steel.

## Delivery program:

CA <sup>®</sup> TIP Pliers		
CA <sup>®</sup> TIP 1	5281.1	1 piece
CA <sup>®</sup> TIP 2	5291.1	1 piece
CA <sup>®</sup> TIP 3	5292.1	1 piece
CA <sup>®</sup> TIP 4	5295.1	1 piece
CA <sup>®</sup> TIP 5	5299.1.	1 piece
CA <sup>®</sup> TIP 6	5323.1	1 piece

Detailed information on use of the pliers can be found here:



SCHEU-DENTAL GmbH  
www.scheu-dental.com

Am Burgberg 20  
58642 Iserlohn · Germany

phone +49 2374 9288-0  
fax +49 2374 9288-90

**SCHEU**  
SCHEU GROUP

