


Stainless steel holder

Product Information



**Important
note:**

To ensure a secure hold of the instruments in the stainless steel holder, tighten the nut if necessary with a suitable pair of pliers (see figure).

Note:

This product contains nickel. In case of nickel allergy make sure to avoid contact with skin/ mucous membrane.

Sterilizable

Mat. Number: 1.4305 / AISI 303

REF 5358

CE



SCHEU-DENTAL GmbH
www.scheu-dental.com

Am Burgberg 20
58642 Iserlohn · Germany

phone +49 2374 9288-0
fax +49 2374 9288-90

SCHEU

SCHEUGROUP

DE/GB 1.000/05/21 G REF PM0253.01

