

CA® Stripping Tools

| REF | Type | Blattstärke mit Beschichtung Blade thickness with coating | Toleranz Tolerance |
|------|-----------------|--|-----------------------|
| 5342 | ● green - left | 0,242 mm | ± 0,02 mm |
| 5343 | ● green - right | 0,217 mm | ± 0,02 mm |
| 5344 | ● green - twin | 0,340 mm | ± 0,02 mm |
| 5345 | ○ grey - left | 0,204 mm | ± 0,02 mm |
| 5346 | ○ grey - right | 0,236 mm | ± 0,02 mm |
| 5347 | ○ grey - twin | 0,276 mm | ± 0,02 mm |
| 5348 | ● red - left | 0,151 mm | ± 0,02 mm |
| 5349 | ● red - right | 0,142 mm | ± 0,02 mm |
| 5350 | ● red - twin | 0,204 mm | ± 0,02 mm |
| 5351 | ○ white - left | 0,097 mm | ± 0,02 mm |
| 5352 | ○ white - right | 0,107 mm | ± 0,02 mm |
| 5353 | ○ white - twin | 0,152 mm | ± 0,02 mm |
| 5354 | ● yellow - twin | 0,118 mm | ± 0,02 mm |

CA® Strip Opener

| REF | Type | Blattstärke mit Beschichtung Blade thickness with coating | Toleranz Tolerance |
|------|-----------------|--|-----------------------|
| 5326 | ● orange - twin | 0,120 mm | ± 0,02 mm |
| 5327 | ● red - twin | 0,170 mm | ± 0,02 mm |

