

<b>REF 3214 - Cutting bur (white)</b> Rough cutting of hard foils	<b>max. 25.000 r/min</b>
<b>REF 3369 - HM-carbide cutter (yellow)</b> Rough cutting of hard foils and hard/soft foils	<b>max. 7.500 r/min</b>
<b>REF 3370 - HM-Tricutter fine (green)</b> Finishing of hard foils and hard/soft foils	<b>max. 7.500 r/min</b>
<b>REF 3377 - Finishing bur (red)</b> Cutting of hard foils and splints with integrated wire elements	<b>max. 7.500 r/min</b>
<b>REF 3226 - Bioperm-Trimmer (blue)</b> Rough trimming of soft foils and hard/soft foils	<b>max. 25.000 r/min</b>
<b>REF 3371 - Polyfix (brown)</b> Finishing of soft foils and hard/soft foils	<b>max. 12.000 r/min</b>
<b>REF 3247 - OSAMU -Polisher (grey)</b> Final polishing of hard/soft foils	<b>max. 15.000 r/min</b>
<b>REF 3376 - DIMO®PRO SLIM (orange)</b> Final trimming and polishing	<b>max. 10.000 r/min</b>

DE/GB 1.000/04/21 G REF 0306.01

phone +49 2374 92 88-0  
fax +49 2374 92 88-90

SCHEU-DENTAL GmbH  
www.scheu-dental.com

**SCHEU**  
SCHEUGROUP