

LC-Primer

Processing Information

Description:

LC- Primer is a light-curing primer for a chemical compound between PET-G pressure moulding material and DURASPLINT® LC and CLEAR-BLOKKER®. It can also be used for the compound of already polymerised DURASPLINT® LC with CLEAR-BLOKKER® in case of build-ups or corrections.

Composition:

Methyl methacrylate, acrylic acide ester, photoinitiator

Curing:

5 minutes in light-curing oven (LC-6 Light oven) equipped with 3 UV-A tubes and 3 VIS tubes (wavelength range 350 – 450 nm), alternating order of tubes.

Packaging unit:

10 ml

Storage condition:

4 – 25°C, not to be exposed to direct sunlight.



SCHEU-DENTAL GmbH
www.scheu-dental.com

phone +49 2374 9288-0
fax +49 2374 9288-90

Am Burgberg 20
58642 Iserlohn · Germany

 Facebook
scheu.dental


Dental Technology

DE/GB 1.500/08/19 G REF PM0084.01
